SAFE | SUSTAINABLE | SMART

SAFER SYRINGE**SAFER STING**

Roncadelle Operations S.R.L.
Via Renolda, 10 - 25030
Castel Mella (BS) - Italy
www.roncadelle-operations.com

ENGLISH**PRODUCT NAME**

Safer® Syringe, safety syringe (hereafter Safer® syringe)
Safer® sting, retractable safety needle (hereafter Safer® needles)

SIZES

Sizes of Syringe: Safety syringe 1ml, 2ml, 2,25ml, 3ml, 5ml.
Needle gauge availability: 20G, 21G, 22G, 23G, 25G, 26G, 27G.
An 18G blunt fill needle is provided separately.

Safer® needle is provided separately. Compatibility between Safer® syringe and Safer® needle is guaranteed for each needle size.

INTENDED USE

Safer® is an automatically Retractable Safety Syringe, it is intended to provide a safe, accurate and reliable method for the aspiration of fluids and for the injection of fluids immediately after filling. It is intended to be used for the administration of medications (intramuscular (IM), intradermal (ID) and subcutaneous (SC)). The Syringe is one assembled with its needle, incorporates a passive safety mechanism which retracts and contains the contaminated needle after use, aiding in the prevention of possible infections due to needle stick injuries and syringe reuse. Needle retraction is activated by the syringe user with one hand.

INSTRUCTION FOR USE

Follow the institution's protocol for applying the aseptic technique (or follow the indication of medical personnel) and skin area preparation. Wash your hands, arrange a suitable, clean place (table) with medicine vial, syringe, alcohol swabs, sharpener and needle.

Peel off the packages to access both the desired Safer® syringe and desired Safer® needles. Attach the Safer® blunt fill needle to the Safer® syringe with a clockwise twist and perform medication preparation and drug aspiration. Then, once the syringe is filled with the desired quantity of medicine, remove the blunt fill needle and dispose of it in a sharps container.

Attach the Safer® Needle for injection to the Safer® syringe with a clockwise twist.

Pull needle cap off the needle and expel air bubbles. Inject medication according to established medical protocol or under the indication of medical personnel. Make sure all of the medication is injected.

Continue to press the plunger button until the retrieval mechanism is activated. The needle will be automatically retracted inside the barrel and a haptic feeling will be perceived as well as a "click" sound.

Do not remove the needle connector after use. Dispose of the syringe in the sharps container.

FEATURES
Biocompatible, sterile, Pyrogen-free.

Harmless to patient and user need once retracted. The needle can be retracted directly from the body.

PRECAUTIONS

* End users should be familiar with instruction for use and be trained in the administration of medication and injections.

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