Language: English

PRODUCT NAME

SafeR® Syringe, safety syringe (SafeR® syringe)

SafeR® Sting, retractable safety needle (SafeR® needles)

SIZES

Sizes of syringe: 1ml, 2ml, 3ml, 5ml.

Needle gauge availability: 20G, 21G, 22G, 23G, 25G, 26G, 27G, 30G.

An 18G blunt fill needle is provided separately.

SafeR® needle is provided separately. Compatibility between SafeR® syringe and SafeR® needle is guaranteed for each needle size.

INTENDED USE

SafeR® is an automatically Retractable Safety Syringe, it is intended to provide a safe, accurate and reliable method for the aspiration of fluids and for the injection of fluids immediately after filling. It is intended to be used for the administration of medications (intramuscular (IM), intradermal (ID) and subcutaneous (SC)). The syringe, once assembled with its needle, incorporates a passive safety mechanism which retracts and contains the contaminated needle after use, aiding in the prevention of possible infections due to needle stick injuries and syringe re-use. Needle retraction is activated by the syringe user with one hand.

INSTRUCTION FOR USE

Follow the institution's protocol for applying the aseptic technique (or follow the indication of medical personnel) and skin area preparation.

Peel off the packages to access both the desired SafeR® syringe and desired SafeR® needle/s. Attach the blunt fill needle to the syringe with a clockwise twist and perform medication preparation and drug aspiration. Then, once the syringe is filled with the desired quantity of medicine, remove the blunt fill needle and dispose of it in a sharps container.

Attach the needle for injection to the syringe with a clockwise twist.

Inject medication, according to established medical protocol or under the indication of medical personnel. Make sure all of the medication is injected.

Continue pressing the plunger button until the retrieval mechanism is activated. The needle will be automatically retracted inside the barrel and a "click" sound will be perceived.

After use, dispose of the device according to the operating procedures of healthcare facilities specific for potentially infectious waste management.

In case of the retrieval mechanism failure, the user may abilitate the mechanism by slightly pulling back the plunger, otherwise dispose of the syringe in the sharps container without re-capping the needle.

Do not remove the needle connector after use.

FEATURES

Biocompatible, Sterile, Pyrogen free.
Safety feature operated with one hand.

PRECAUTIONS

- End users should be familiar with instruction for use and be trained in the administration of medication and injections.

Document Owner	Michael Rota Scalabrini
----------------	-------------------------

- Do not use if single unit package is damaged or needle cap is detached.
- The syringes and needles are sterile (EtO sterilization) until removed from their protective packaging. They are single use only.
- The needle retrieval mechanism can only be activated when the medication or fluid is fully injected.
- Do not push the syringe further inside the body while activating the retracting mechanism.
- Discard of the syringe immediately after use without any further utilization.

WARNINGS

- SafeR® retractable needle is designed for use with SafeR® syringes only.
- Syringe re-use could lead to infection or transmission of diseases: do not try to re-use the device
- Device shall not be used on patient having a body weight lower than 3 kg.
- Do not use this product to draw blood or for endovenous punctures.
- SafeR® is a single use device: do not re-sterilize.
- Blunt fill needle is designed only for drug preparation and it is not intended to come in contact with the human body. DO NOT USE IT FOR MEDICATION INJECTION
- Report any device malfunction, deterioration in device performance, inadequate instructions, or inadequate labeling to the manufacturer and competent authority, as if not reported, this may result in death or serious injury, or may lead to death or serious deterioration in state of health.

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